DÉSCRIPTION / DENOMINATION / سُبَقَة	RÉFÉRENCE / REFERENCE / المرجع
Stéthoscope MAGISTER II Adulte simple pavillon Silver MAGISTER II Il adulto simple estetoscopio - Silver Estetoscopio MAGISTER II para adultos simple campana plata MAGISTER II Erwachsenen-Einzellokopfstethoskop - Silber Stetoscopia a padiglione singolo MAGISTER II - Argento MAGISTER-II-stethoscoop voor volwassenen met één kop - zilver MAGISTER II akuskin yhdys päänt stetoskope - Hopea Täikasavutte ue peaga stetoskoop MAGISTER II - Hõbe MAGISTER II saugusujų vienos galvos nėrauskopas - Sidabrinė	507 250
Stéthoscope MAGISTER II Adulte double pavillon Black MAGISTER II Adult stethoscope dual chestpiece - Black Estetoscopio MAGISTER II para adultos doble campana negro MAGISTER II Erwachsenen-Einzellokopfstethoskop - Schwarze Stetoscopia a padiglione doppio MAGISTER II - Nero MAGISTER-II-stethoscoop voor volwassenen met één kop - zwart MAGISTER II akuskin yhdys päin stetoskope - Musta Täikasavutte ue peaga stetoskoop MAGISTER II - Must MAGISTER II saugusujų dvigubu vėliau - Juodas	507 260
Stéthoscope MAGISTER II Adulte double pavillon Silver MAGISTER II Adult stethoscope dual chestpiece - Silver Estetoscopio MAGISTER II para adultos doble campana plata MAGISTER II Erwachsenen-Einzellokopfstethoskop - Silber Stetoscopia doppio MAGISTER II per adulti - Argento MAGISTER-II-stethoscoop Volwassen dubbele vlag - Zilver MAGISTER-II-stethoscoop alkuisten kaksiosainen topeltipäät - Hopea MAGISTER-II-saugusujų dvigubu vėliau - Sidabrinė	507 251
Stéthoscope MAGISTER II Adulte double pavillon Black MAGISTER II Adult stethoscope dual chestpiece - Black Estetoscopio MAGISTER II para adultos doble campana negro MAGISTER II Erwachsenen-Einzellokopfstethoskop - Schwarze Stetoscopia a padiglione doppio MAGISTER II - Nero MAGISTER-II-stethoscoop voor volwassenen met één kop - zwart MAGISTER II akuskin yhdys päin stetoskope - Musta Täikasavutte ue peaga stetoskoop MAGISTER II - Must MAGISTER II saugusujų dvigubu vėliau - Juodas	507 261

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USER'S MANUAL		
Thank you for choosing the MAGISTER II stethoscope.		
Choosing Spenger's instrument is trusting the experience and recognized excellence reconduced since 1907.		
MAGISTER II is the resulting from many years of research and development aiming at providing the best quality, design and comfort to health professionals.		
This medical device is exclusively manufactured in accordance with the Directive 93/42 EEC, relating to medical devices and is subject to quality controls.		
INSTRUCTIONS FOR USE		
Please read carefully the user manual before using the product and keep it safe for further reference.		
INTENDED MEDICAL USE		
MAGISTER II is a professional equipment (can be used by trained health professionals) for medical diagnostic. It may be used for heart auscultation, lungs... as well as measuring the blood pressure.		
Caution: This device must be used according to the instructions		
INSTRUCTIONS FOR USE		
Fitting and adjusting the binaural tubes		
The binaural tubes are designed to automatically adjust to your ear's listening angle. Ensure that the branches of the binaural tubes are pointing forwards, thus matching the direction of the auditory canal [fig. 1].		
The tension of the spring can be adjusted for improved comfort and better adjustment to your morphology.		
To increase the spring tension, grasp the branches of the binaural tubes at the level of the earpieces and cross them over [fig. 2].		
To reduce the spring tension, grasp the branches at the base of the binaurals and gradually move them apart to complete extension [fig. 3].		
Selecting chestpiece for MAGISTER II DUAL		
Thanks to the innovative chest-piece indicator, you can see which side of the dual head is activated.		
Using the multi-frequency floating diaphragm		
The MAGISTER II stethoscope enables you to perform various types of auscultation. The multi-frequency floating diaphragm can be used to listen to both low and high frequency sounds without the need to flip the chestpiece.		
• Low frequencies:		
Apply very mild pressure to the chestpiece. The diaphragm is held in place by a PVC ring, enabling low-frequency sounds to resonate.		
• High frequencies:		
Apply firm pressure on the chestpiece. The diaphragm makes contact with an internal ring, thus allowing high-frequency sounds to be heard.		
To alternate between listening to high and low-frequency sounds, simply alter the pressure applied to the chestpiece.		
Using the small bell chestpiece		
The small bell chestpiece allows listening to low frequencies sounds.		
MAINTENANCE		

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